


CLINICON LAVAGE TRAY

Inflatable surgical tray designed to prevent spillage of blood and lavage fluids during surgery. Satisfies COSHH regulations 2002 to reduce the risks of cross-infection of patients and staff during transfer to and within the hospital environment.



Simple and effective.

Cushioned walls to protect patients and staff.

Inflated using standard air line and universal tubing.

Drainage port for connection to universal tubing and standard vacuum suction unit.

Sieve device to ensure rapid drainage of infectious tissue fluids.

Self-adhesive tapes allow the tray to be folded for safe disposal after use.

Reduced need for additional staff within the operating field.

Cost effective for rapid patient turnaround times.

Patent Protected.

Supplied sterile, in dispensing packs of 10 units.



Code	Description	Dimensions
CL 460-800	Clinicon Lavage tray	80 x 45cm (flat)