

### LEFT SIDED SET

Endobronchial tube, connector set and two control valve suction catheters with graduation mark. Single use. Supplied sterile.

- Medical grade plastic stylet
- Two soft low-pressure cuffs
- Atraumatic rounded tip
- X-ray detectable line

Size	Code	NPC Code
26ch	PRS-26L	-
28ch	PRS-28L	<b>DFD1131</b>
32ch	PRS-32L	<b>DFD1132</b>
35ch	PRS-35L	<b>DFD1133</b>
37ch	PRS-37L	<b>DFD1134</b>
39ch	PRS-39L	<b>DFD1135</b>
41ch	PRS-41L	<b>DFD1136</b>



AVAILABLE via  
NHS SUPPLY CHAIN

### RIGHT SIDED SET

Endobronchial tube, connector set and two control valve suction catheters with graduation mark. Single use. Supplied sterile.

- Medical grade plastic stylet
- Two soft low-pressure cuffs
- Atraumatic rounded tip
- Murphy eye
- X-ray detectable line

Size	Code	NPC Code
26ch	PRS-26R	-
28ch	PRS-28R	<b>DFD1137</b>
32ch	PRS-32R	<b>DFD1138</b>
35ch	PRS-35R	<b>DFD1139</b>
37ch	PRS-37R	<b>DFD1140</b>
39ch	PRS-39R	<b>DFD1141</b>
41ch	PRS-41R	<b>DFD1142</b>

