

TRACHEAL TUBE WITHOUT CUFF

Transparent thermosensitive PVC.
Oral/Nasal. Murphy type, with Xray line.

Sterile, single use. Supplied in boxes of 10.

ID(mm)	OD(mm)	Code (PVC)	Code (Siliconised)
2.0	3.0	22-2000	22-2010
2.5	3.5	22-2500	22-2510
3.0	4.2	22-3000	22-3010
3.5	4.8	22-3500	22-3510
4.0	5.4	22-4000	22-4010
4.5	6.0	22-4500	22-4510
5.0	6.9	22-5000	22-5010
5.5	7.5	22-5500	22-5510
6.0	8.2	22-6000	22-6010
6.5	8.8	22-6500	22-6510
7.0	9.7	22-7000	22-7010
7.5	10.3	22-7500	22-7510
8.0	11.0	22-8000	22-8010
8.5	11.7	22-8500	22-8510
9.0	12.3	22-9000	22-9010
9.5	13.0	22-9500	22-9510
10.0	13.7	22-1000	22-1010



TRACHEAL TUBE WITHOUT CUFF, SOFT

Oral/Nasal. Good for long-term ventilation. Soft flexible plastic reduces risk of trauma. Double graduation mark every cm. Murphy type, with Xray line.

Sterile, single use. Supplied in boxes of 10.

ID(mm)	OD(mm)	Code (PVC)	Code (Siliconised)
2.0	3.3	23-2000	23-2010
2.5	4.0	23-2500	23-2510
3.0	4.5	23-3000	23-3010
3.5	5.1	23-3500	23-3510
4.0	5.8	23-4000	23-4010
4.5	6.5	23-4500	23-4510
5.0	7.3	23-5000	23-5010
5.5	8.0	23-5500	23-5510
6.0	8.6	23-6000	23-6010
6.5	9.2	23-6500	23-6510
7.0	9.7	23-7000	23-7010
7.5	10.3	23-7500	23-7510
8.0	11.0	23-8000	23-8010
8.5	11.7	23-8500	23-8510
9.0	12.3	23-9000	23-9010
9.5	13.0	23-9500	23-9510
10.0	13.7	23-1000	23-1010



TRACHEAL TUBE WITH LOW PRESSURE CUFF

Transparent thermosensitive PVC.
Oral/Nasal. Available Magill or Murphy type, with Xray line.

Sterile, single use. Supplied in boxes of 10.

ID(mm)	OD(mm)	Cuff(mm)	Code (Magill type)		Code (Murphy type)	
			PVC		PVC	Siliconised
2.5	3.5	8	32-2500		35-2500	35-2510
3.0	4.2	8	32-3000		35-3000	35-3010
3.5	4.8	10	32-3500		35-3500	35-3510
4.0	5.4	10	32-4000		35-4000	35-4010
4.5	6.0	12	32-4500		35-4500	35-4510
5.0	6.9	14	32-5000		35-5000	35-5010
5.5	7.5	17	32-5500		35-5500	35-5510
6.0	8.2	20	32-6000		35-6000	35-6010
6.5	8.8	20	32-6500		35-6500	35-6510
7.0	9.7	23	32-7000		35-7000	35-7010
7.5	10.3	25	32-7500		35-7500	35-7510
8.0	11.0	26	32-8000		35-8000	35-8010
8.5	11.7	27	32-8500		35-8500	35-8510
9.0	12.3	29	32-9000		35-9000	35-9010
9.5	13.0	29	32-9500		35-9500	35-9510
10.0	13.7	32	32-1000		35-1000	35-1010



TRACHEAL TUBE WITH HIGH VOLUME CUFF

Transparent siliconised PVC. Oral/Nasal. Thin-walled high volume cuff and marked pilot balloon. Good for long-term ventilation. Murphy eye, with Xray line.

Sterile, single use. Supplied in boxes of 10.

ID(mm)	OD(mm)	Cuff(mm)	Code
5.0	6.9	20	24-5095
5.5	7.5	20	24-5595
6.0	8.2	25	24-6095
6.5	8.8	25	24-6595
7.0	9.7	28	24-7095
7.5	10.3	28	24-7595
8.0	11.0	30	24-8095
8.5	11.7	30	24-8595
9.0	12.3	34	24-9095
9.5	13.0	34	24-9595
10.0	13.7	34	24-1095





REINFORCED TRACHEAL TUBE

Oral/Nasal. Soft flexible plastic with stainless steel kink-resistant reinforcement. Available with or without cuff. Murphy type, with marked pilot balloon.

Sterile, single use. Supplied in boxes of 10.

ID(mm)	OD(mm)	Cuff(mm)	Code (with cuff)		Code (without cuff)	
			PVC	Siliconised	PVC	Siliconised
3.0	4.8	8	33-3000	33-3010	26-3000	26-3010
3.5	5.4	10	33-3500	33-3510	26-3500	26-3510
4.0	6.0	10	33-4000	33-4010	26-4000	26-4010
4.5	6.9	12	33-4500	33-4510	26-4500	26-4510
5.0	7.5	14	33-5000	33-5010	26-5000	26-5010
5.5	8.2	17	33-5500	33-5510	26-5500	26-5510
6.0	8.8	20	33-6000	33-6010	26-6000	26-6010
6.5	9.6	20	33-6500	33-6510	26-6500	26-6510
7.0	10.2	23	33-7000	33-7010	26-7000	26-7010
7.5	11.0	25	33-7500	33-7510	26-7500	26-7510
8.0	11.6	26	33-8000	33-8010	26-8000	26-8010
8.5	12.2	27	33-8500	33-8510	26-8500	26-8510
9.0	12.9	29	33-9000	33-9010	26-9000	26-9010
9.5	13.6	29	33-9500	33-9510	26-9500	26-9510
10.0	14.2	32	33-1000	33-1010	26-1000	26-1010



TRACHEAL TUBE WITH SUCTION LUMEN

Siliconised PVC. Oral/Nasal. Allows suction of secretions from above the cuff, minimising risk of infection. Delicate low-pressure cuff. Suction lumen with cap or connector. Murphy eye, with Xray line and marked pilot balloon.

Sterile, single use. Supplied in boxes of 10.

ID(mm)	OD(mm)	Cuff(mm)	Code (with cap)	Code (with connector)
5.0	6.9	14	11-5015	12-5011
5.5	7.5	17	11-5515	12-5511
6.0	8.2	20	11-6015	12-6011
6.5	8.8	20	11-6515	12-6511
7.0	9.7	23	11-7015	12-7011
7.5	10.3	25	11-7515	12-7511
8.0	11.0	26	11-8015	12-8011
8.5	11.7	27	11-8515	12-8511
9.0	12.3	29	11-9015	12-9011
9.5	13.0	29	11-9515	12-9511
10.0	13.7	32	11-1015	12-1011



MICROLARYNGEAL TRACHEAL TUBE

Thin-walled, high volume, low-pressure cuff, for use in laryngeal surgery. Oral/Nasal. Murphy type, with Xray line.

Sterile, single use. Supplied in boxes of 10.

ID(mm)	OD(mm)	Cuff(mm)	Length(mm)	Code (PVC)	Code (Siliconised)
4.0	5.4	20	368	10-4000	10-4010
5.0	6.9	25	368	10-5000	10-5010
6.0	8.2	28	368	10-6000	10-6010

PREFORMED TRACHEAL TUBE, NASAL, NORTH

Soft flexible PVC available clear or opaque and with or without low-pressure cuff. Murphy type, with Xray line and marked pilot balloon.

Sterile, single use. Supplied in boxes of 10.


ID(mm)	OD(mm)	Cuff(mm)	CLEAR				OPAQUE			
			Code (with cuff)		Code (without cuff)		Code (with cuff)		Code (without cuff)	
			PVC	Siliconised	PVC	Siliconised	PVC	Siliconised	PVC	Siliconised
3.0	4.2	8	38-3000	38-3010	40-3000	40-3010	38-3004	38-3014	41-3004	41-3014
3.5	4.8	10	38-3500	38-3510	40-3500	40-3510	38-3504	38-3514	41-3504	41-3514
4.0	5.4	10	38-4000	38-4010	40-4000	40-4010	38-4004	38-4014	41-4004	41-4014
4.5	6.0	12	38-4500	38-4510	40-4500	40-4510	38-4504	38-4514	41-4504	41-4514
5.0	6.9	14	38-5000	38-5010	40-5000	40-5010	38-5004	38-5014	41-5004	41-5014
5.5	7.5	17	38-5500	38-5510	40-5500	40-5510	38-5504	38-5514	41-5504	41-5514
6.0	8.2	20	38-6000	38-6010	40-6000	40-6010	38-6004	38-6014	41-6004	41-6014
6.5	8.8	20	38-6500	38-6510	40-6500	40-6510	38-6504	38-6514	41-6504	41-6514
7.0	9.7	23	38-7000	38-7010	40-7000	40-7010	38-7004	38-7014	41-7004	41-7014
7.5	10.3	25	38-7500	38-7510	40-7500	40-7510	38-7504	38-7514	41-7504	41-7514
8.0	11.0	26	38-8000	38-8010	40-8000	40-8010	38-8004	38-8014	41-8004	41-8014

PREFORMED TRACHEAL TUBE, ORAL, SOUTH

Available with or without low-pressure cuff. Murphy type, with Xray line and marked pilot balloon.

Sterile, single use. Supplied in boxes of 10.

ID(mm)	OD(mm)	Cuff(mm)	Code (with cuff)		Code (without cuff)	
			PVC	Siliconised	PVC	Siliconised
3.0	4.2	8	34-3000	34-3010	37-3000	37-3010
3.5	4.8	10	34-3500	34-3510	37-3500	37-3510
4.0	5.4	10	34-4000	34-4010	37-4000	37-4010
4.5	6.0	12	34-4500	34-4510	37-4500	37-4510
5.0	6.9	14	34-5000	34-5010	37-5000	37-5010
5.5	7.5	17	34-5500	34-5510	37-5500	37-5510
6.0	8.2	20	34-6000	34-6010	37-6000	37-6010
6.5	8.8	20	34-6500	34-6510	37-6500	37-6510
7.0	9.7	23	34-7000	34-7010	37-7000	37-7010
7.5	10.3	25	34-7500	34-7510	37-7500	37-7510
8.0	11.0	26	34-8000	34-8010	37-8000	37-8010
8.5	11.7	27	34-8500	34-8510	37-8500	37-8510
9.0	12.3	29	34-9000	34-9010	37-9000	37-9010



PREFORMED TRACHEAL TUBE, ORAL, NORTH

Available with or without low-pressure cuff. Murphy type, with Xray line and marked pilot balloon.

Sterile, single use. Supplied in boxes of 10.

ID(mm)	OD(mm)	Cuff(mm)	Code (with cuff)		Code (without cuff)	
			PVC	Siliconised	PVC	Siliconised
3.0	4.2	8			37-3027	37-3037
3.5	4.8	10			37-3527	37-3537
4.0	5.4	10			37-4027	37-4037
4.5	6.0	12			37-4527	37-4537
5.0	6.9	14	34-5027	34-5037	37-5027	37-5037
5.5	7.5	17	34-5527	34-5537	37-5527	37-5537
6.0	8.2	20	34-6027	34-6037	37-6027	37-6037
6.5	8.8	20	34-6527	34-6537	37-6527	37-6537
7.0	9.7	23	34-7027	34-7037	37-7027	37-7037
7.5	10.3	25	34-7527	34-7537	37-7527	37-7537
8.0	11.0	26	34-8027	34-8037	37-8027	37-8037
8.5	11.7	27	34-8527	34-8537	37-8527	37-8537
9.0	12.3	29	34-9027	34-9037	37-9027	37-9037

