

LARYNGEAL MASKS, ECONOMY

- Soft PVC with atraumatic translucent silicone cuff.
- MRI compatible. Contains Phthalates.

Single use. Supplied sterile in boxes of 10.

Size	Suitable	Code
1	up to 5kg	501000
1.5	5 - 10kg	501500
2	10 - 20kg	502000
2.5	20 - 30kg	502500
3	30 - 50kg	503000
4	50 - 70kg	504000
5	70 - 100kg	505000
6	over 100kg	506000



LARYNGEAL MASKS, ECONOMY, ALL PVC

- All PVC including cuff.
- Non MRI compatible. Contains Phthalates.

Single use. Supplied sterile in boxes of 10.

Size	Suitable	Code
1	up to 5kg	551000
1.5	5 - 10kg	551500
2	10 - 20kg	552000
2.5	20 - 30kg	552500
3	30 - 50kg	553000
4	50 - 70kg	554000
5	70 - 100kg	555000
6	over 100kg	556000



LARYNGEAL MASKS, SILICONE

- All silicone with atraumatic silicone cuff.
- MRI compatible.

Available single use (translucent cuff) or reusable (pink cuff).
Supplied sterile in boxes of 10.

Size	Suitable	Code (single use)	Code (reusable)
1	up to 5kg	511000	521000
1.5	5 - 10kg	511500	521500
2	10 - 20kg	512000	522000
2.5	20 - 30kg	512500	522500
3	30 - 50kg	513000	523000
4	50 - 70kg	514000	524000
5	70 - 100kg	515000	525000
6	over 100kg	516000	526000



LARYNGEAL MASKS, FLEXIBLE

- Kink-resistant flexible silicone tube with atraumatic silicone cuff.
- Non MRI compatible.

Available single use (translucent cuff) or reusable (pink cuff).
Supplied sterile in boxes of 10.

Size	Suitable	Code (single use)	Code (reusable)
1	up to 5kg	531000	541000
1.5	5 - 10kg	531500	541500
2	10 - 20kg	532000	542000
2.5	20 - 30kg	532500	542500
3	30 - 50kg	533000	543000
4	50 - 70kg	534000	544000
5	70 - 100kg	535000	545000
6	over 100kg	536000	546000

