


- Original Dr Frank Robertshaw design.
- “D” shaped lumina provides maximum cross sectional area, minimising resistance to gas flows.
- Blue bronchial cuff assists safe bronchoscopic positioning.
- Disposable connection set included.
- Environmentally friendly, natural biodegradable materials.

Supplied sterile, individually packed, with connector set.

LEFT SIDED TUBE

Sterile, single use.

Size	Code
Extra small	SD1200
Small	SD1201
Medium	SD1202
Large	SD1203



RIGHT SIDED TUBE

Sterile, single use.

Size	Code
Extra small	SD1300
Small	SD1311
Medium	SD1312
Large	SD1313

