



**P3 Medical Limited**  
**1 Newbridge Close, Bristol BS4 4AX**  
**United Kingdom**

Phone +44 (0)117 972 8888  
 Fax +44 (0)117 972 4863  
 Email: info@p3-medical.com

**Nelaton Catheter**

General description

The Nelaton catheter is a fine hollow red rubber re-usable tube which is inserted into the bladder to drain urine away.


Intended Use

Inserted into the bladder to drain urine away.

Technical details

Materials: Natural Red Rubber

Latex status: Contains Latex

Description	Size (Ch)	Part Number
<b>Nelaton Catheter Non Sterile Reusable</b>  	8Ch	700108PHX *
	10Ch	700110PHX *
	12Ch	700112PHX *
	14Ch	700114PHX *
	16Ch	700116PHX *
	18Ch	700118PHX *
	20Ch	700120PHX *
	22Ch	700122PHX *
	24Ch	700124PHX *
	26Ch	700126PHX *
	28Ch	700128PHX *
	30Ch	700130PHX *

\*Suffix code PHX may vary for own brand products

Packaging

Packed in 50's or 20's dependent on size and sealed poly bag to be sold non sterile.

Reprocessing Limitations

Frequent (thermal) reprocessing affects the service life of these products (up to 24 reprocessing cycles). The end of the product service life can also be limited by wear and damage due to use

Cleaning

Nelaton Catheters may be cleaned with mild detergents e.g. Tri-Sodium Phosphate or a neutral liquid soap and then dried at ambient temperature

Regulatory compliance

Meets the essential requirements of directive 93/42/EEC as amended 2007/47/EC for Class I medical devices

[www.p3-medical.co.uk](http://www.p3-medical.co.uk)

Registered number 1072913 England  
 VAT registration number GB 138 5525 56