

TRACHESEAL

Post-decannulation tracheostomy wound dressing, for faster healing and greater patient comfort following removal of tracheostomy tube.



- Unique dressing specifically designed to address stoma closure difficulties.
- Centre button plate enables accurate and consistent location of the tracheostomy wound.
- Gentle pressing of the button plate prevents unwanted passage of air through the wound during actions such as talking, coughing and sneezing.
- Proven hydrocolloid dressing material absorbs exudate and has bacteriostatic action.
- Promotes faster healing and restoration of speech.
- Improves whole patient experience following tracheostomy tube decannulation, leading to shorter hospital stays and fewer complications.

Individually packed sterile, 10 units per carton.

Code	Description	Dimensions
TS0001	Tracheostomy wound dressing	100mm dia

Patents pending